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

Organization	Ultrapure water quality/grade
European Pharmacopoeia	Purified water in bulk
U.S. Pharmacopoeia	Purified water in bulk
Japanese Pharmacopoeia	Purified water
Chinese Pharmacopoeia	Purified water
ASTM® D1193	Type I water
ISO 3696	Grade 1 water
Chinese National Standard GB/T 6682	Grade 1 water
Chinese National Standard GB/T 33087	Ultrapure water
JIS K 0557	A4 water
Clinical and Laboratory Standards Institute (CLSI®)	Clinical Laboratory Reagent Water (CLRW)

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Parameter	Milli-Q® TOC Indicator
Monitoring frequency	At dispense
Accuracy	Accurate indication within the range
TOC values display	≤ 5 ppb, if 0–5 ppb ≤ 10 ppb, if 6–10 ppb > 10–999 ppb, a whole number is displayed
TOC measuring process	Inline, post-dispense
Photooxidation UV lamp	Low pressure mercury lamp, 185 nm
UV lamp replacement frequency	Once every 2 years

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







Q-POD® dispenser

Tubing and port requirements

Parameter	Description
Feed water connection	½" GAZ
Distance from feed water port	Maximum 5 m
Distance from purification unit to wall-mounted Q-POD® dispenser	Maximum 3 m
Distance from purification unit to tank	Maximum 5 m
Power entry	Connection IEC13
ON/OFF switch	Available on the unit
Water sensor port	Fits with standard leak detector (Cat. No. ZWATSENA1)
Ethernet port	IEEE P802.3

Electrical connections and specifications

Parameter	Description
Power source voltage	100~240 VAC ± 10%
Power frequency	50/60 Hz ± 2 Hz
Power used	Maximum 200 VA
Power cord length	2.5 m, plug: IEC13 female
Operational temperature	4~40°C
Altitude	Up to 3000 m
Tank material	High purity polyethylene

Weights

System type		Dry weight	Shipping weight	Operating weight
Purification unit + HMI touchscreen	7008	19.0 kg*	22.4 kg	25.2 kg*
	7016	19.7 kg*	23.1 kg	26.5 kg*
HMI touchscreen		0.58 kg	**	0.58 kg
HMI touchscreen mounting kits	Benchtop	0.42 kg	0.56 kg	0.42 kg
	Wall	0.49 kg	1.4 kg	0.49 kg
Q-POD® dispenser	System	1.2 kg	1.5 kg	1.2 kg
	Wall	2.2 kg	2.8 kg	2.2 kg
Storage Tank	25 L	6.7 kg	8.5 kg	31.7 kg
	50 L	7.6 kg	10.6 kg	57.6 kg
	100 L	10.9 kg	12.8 kg	110.9 kg

* HMI touchscreen not included. ** Included in shipping box with purification unit.

HMI touchscreen specifications

Screen descriptions and functionalities

Parameter	Description
Capacitive touchscreen	Size: 7" (18 cm); Resolution: 800 × 480
USB port	USB 2.0 Highspeed standard
Display in 9 languages	Chinese / English / French / German / Italian / Japanese / Portuguese / Russian / Spanish

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Feed water requirements

Feed water quality	Potable tap water
Pressure	1~6 bar
Temperature	5~35°C
Conductivity	< 2000 μ S/cm at 25°C
Dissolved CO ₂	< 30 ppm
Free chlorine	< 3 ppm*
Fouling Index (FI)	< 10**
pH	4~10
Total Organic Carbon (TOC)	< 2 ppm
Langlier Saturation Index (LSI)	< 0.3
Hardness (as CaCO ₃)	< 300 ppm
Silica	< 30 ppm

* If free chlorine between 1 and 3 ppm, apply PrePak PRPK00001. ** If FI between 5 and 10, apply PrePak PRPK000A1.

Ultrapure, Type 1 water specifications¹ (from Q-POD® dispenser)

Resistivity ²	18.2 M Ω · cm at 25°C
Conductivity	0.055 μ S/cm at 25°C
TOC	≤ 5 ppb
Particles ³	No particles with size > 0.22 μ m
Bacteria	< 0.01 cfu/mL (< 10 cfu/L) ⁴ < 0.005 cfu/mL (< 5 cfu/L) ⁵
Pyrogens (endotoxins) ⁶	< 0.001 EU/mL
RNases ⁶	< 1 pg/mL
DNases ⁶	< 5 pg/mL
Proteases ⁶	< 0.15 μ g/mL
Flow rate	< 2 L/min

1 These values are typical and may vary depending on the nature and concentration of contaminants in the feed water.

2 Resistivity can also be displayed non-temperature-compensated as required by USP.

3 With Millipak® or Millipak® Gold filter.

4 With Millipak® or Biopak® filter.

5 With Millipak® Gold filter when installed and used in a laminar flow hood.

6 With Biopak® polisher.

Tank water / Type 3 water delivery

Dispensing tank water is possible provided that a tank front valve is installed.

Reverse osmosis (RO), Type 3 water specifications

Resistivity	> 0.05 M Ω · cm at 25°C
RO ionic rejection	97~98%
Organics rejection	≥ 99% (depending on the type of molecule)
TOC	< 200 ppb
Colloids	< 1000 ppb
Bacteria	< 1000 cfu/mL (with ASM option installed)
Production flow rate	8 L/h (Milli-Q® EQ 7008) 16 L/h (Milli-Q® EQ 7016)

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Water purification systems, dispensers and mounting kits	Catalog number
Milli-Q® EQ 7008 system (8 L/h production flow rate) & HMI touchscreen	ZEQ7008T0C
Milli-Q® EQ 7016 system (16 L/h production flow rate) & HMI touchscreen	ZEQ7016T0C
Benchtop kit for HMI touchscreen	BTEQ0DKT
Wall mounting kit for HMI touchscreen	WMEQ0DKT
System-mounting kit for Q-POD® dispenser (includes Q-POD® unit)	SMEQ00KT
Wall-mounting kit for Q-POD® dispenser (includes Q-POD® unit)	WMEQ0RKT

Water storage tanks	Catalog number
Milli-Q® storage tank, 25 L	TANKA025
Milli-Q® storage tank, 50 L	TANKA050
Milli-Q® storage tank, 100 L	TANKA100
Milli-Q® storage tank top assembly (includes ASM)	TANKT0PA1
Milli-Q® storage tank top assembly (no ASM)	TANKT0PEQ

Purification consumables	Catalog number
IPAK Gard®/IPAK Meta®/IPAK Quanta®/vent filter consumable kit	EQ70XXPKT1
IPAK Gard®/IPAK Meta®/IPAK Quanta®/vent filter consumable kit for hard water	EQ70XXPKT1H
IPAK Meta®/IPAK Quanta® consumable kit	IPAKKITA1
IPAK Gard® pretreatment pack	IPAKGARA1
IPAK Gard® pretreatment pack for hard water	IPAKGARH1
Vent filter	TANKV01A1
Vent filter HF (for high-flow applications)	TANKVH1A1
UV lamp	ZEQ7UVLP0

Application POD-Paks	Catalog number
Millipak® 0.22 µm filter	MPGP002A1
Millipak® Gold 0.22 µm sterile filter	MPGPG02A1
Biopak® polisher	CDUFBI0A1
LC-Pak® polisher	LCPAK00A1
EDS-Pak® polisher	EDSPAK0A1
VOC-Pak® polisher	V0CPAK0A1

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Accessories & Connectors	Catalog number
System wall mounting bracket	SYSTFIXA1
Tank wall mounting bracket	TANKFIXA1
Tank valve kit	ZFTVK07A1
Connector 2 m system-to-storage tank	ZFC0NN2ST
Connector 5 m system-to-storage tank	ZFC0NN5ST
Water sensor	ZWATSENA1
Foot pedal	ZMQSFTSA1
Alarm relay cable	ZMQ0ALCA1
External solenoid valve for feed water	EXTSV00A1
Washer distribution kit 230 V (right)	ZWDK5R100
Washer distribution kit 230 V (left)	ZWDK5L100
Washer distribution kit 115 V (right)	ZWDK6R100
Washer distribution kit 115 V (left)	ZWDK6L100
Washer distribution kit adaptor	ZWDKADPA1
Wall mounting bracket for washer distribution kit	WMBWASH1
Multi system installation kit	ZIQ7MSKT1
System care	Catalog number
ROCare A - Acidic care	ZWACID012
ROCare B - Basic care	ZWBASE012
ROProtect C - Chlorine tablets	ZWCL01F50
EfferSan Effervescent Tablets (USA)	5874316024
EfferSan Effervescent Tablets (CAN)	5874316024C

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© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Milli-Q, MyMilli-Q, Q-POD, Elix, IPAK Gard, IPAK Meta, IPAK Quanta, IQnano, Jetpore, Millipak, Biopak, VOC-Pak, EDS-Pak and LC-Pak은 Merck KGaA, Darmstadt, Germany 또는 그 계열사의 상표입니다. 기타 모든 상표는 해당 소유자의 자산입니다. 상표에 관한 자세한 정보는 공개적으로 열람 가능한 출처를 통해 이용할 수 있습니다.