

Cell Strainer 런칭기념 이벤트



40µm / 70µm / 100µm Fits into
50ml Centrifuge Tubes



Gamma radiation sterilized

Features:

- Made of sturdy polypropylene and strong nylon mesh
- Three different pore sizes available: 40µm, 70µm, or 100µm, strainers are color-coded according to pore size for convenience
- Fits perfectly into 50ml Centrifuge Tubes
- Gamma radiation sterilized, non-pyrogenic



Cat #	Size	Color	Sterile	Unit	행사가
CBM-CARE151040	40µm	Purple	Yes	100/case	100,000
CBM-CARE151070	70µm	Orange	Yes	100/case	100,000
CBM-CARE151100	100µm	Green	Yes	100/case	100,000

GOLD PACIFIC
골드퍼시픽

(주)골드퍼시픽 강동지점
서울시 강동구 올림픽로 834 한강씨티라이프빌딩 2층
문의 ☎ 02-427-3669 ✉ curebio@curebio.co.kr