


실험실의

SUPERHERO

최강의 어벤져스들이 모였다!!

할인 혜택을 누리보세요



할인 구매와 더불어 추가적인
혜택을 받아보실 수 있습니다



50만원 이상 구매시
샤오미 미밴드

행사기간 : 6월 7일 ~ 6월 30일


Gibco FBS Fingerprinting

어떤 FBS를 신뢰하시나요?

당사는 혈청을 공급하는 최고의 글로벌 기업으로서 고객이 Gibco와 같은 공급업체에게 요구하고 기대하는 제품과 최고 품질의 혈청을 확보하기 위해 다각적인 노력을 펼치고 있음을 고객에게 알려 신뢰도를 높이고자 합니다.

Gibco FBS Fingerprinting 서비스로 FBS 원산지를 최대한 보장해 드립니다.

뉴질랜드에 본사를 둔 서비스 제공업체 Ortain은 지문 검사를 실시합니다. 이 방법은 육류, 꿀, 우유, 치즈, 와인 제품의 원산지 증명과 위조 방지 용도로 성공적으로 사용되어 왔습니다. 원소는 환경, 토양, 식물, 물에서 자연적으로 발생하고, 그 수치와 유형은 지역에 따라 다양합니다. 동물은 환경에서 이 자연적인 화학적 특성을 흡수하여 원산지와 연관되는 고유한 지화학적 프로필을 나타냅니다. Ortain은 혈청에서 이 원소를 분석함으로써 당사 FBS의 원산지를 확인할 수 있습니다.



Gibco FBS는 업계 최고 수준의 보증인 화학적 Fingerprinting 서비스를 제공합니다. Fingerprinting 서류를 신청하세요!



Gibco™ FBS Exosome-depleted FBS

Bovine exosome을 제거하기위한 FBS 초원심분리 과정을 18시간씩 소요하셨나요? 그럼에도 불구하고 일반적으로 50%정도로만 Exosome이 제거될 뿐이고, 또한 FBS안의 세포 성장에 필요한 일부 영양소가 감소되어 기존과 일치하지 않는 세포성장을 야기시키는 결과를 받아보셨나요?

여기 Exosome을 90%이상 제거하고, 세포성장에도 일관성있으며 고성능을 발휘하는 Gibco Exosome-depleted FBS를 소개합니다.

지금 바로 문의하세요~ 특별가로 드립니다.

Gibco ExpiCHO Expression System

Transient 단백질 발현의 미래로 초대합니다. Transient 시스템으로 최대 3g/L 생산



ExpiCHO Expression System에는 발현율이 높은 CHO cell line, 화학적으로 규명된 (chemically defined) 동물 유래 성분이 없는 배양 배지, 최적화된 배양 Feeds, 고효율 transfection reagent 등이 모두 들어있어 상승 작용을 통해 Gibco™ FreeStyle™ CHO Expression System보다 무려 160배나 더 높은 효율과 Gibco™ Expi293™ Expression System보다 4배 높은 효율을 확보할 수 있습니다. human IgG 단백질에서 최대 3 g/L의 발현 수준을 자랑합니다.

이제 연구 초기 단계부터 CHO cell로 시작하세요.

Cell Therapy Systems (CTS)

줄기세포의 기초연구부터 임상연구까지 지원하는 유용한 제품



줄기세포 치료 연구를 처음부터 제대로 시작하기 위해서는, 임상 적용을 위한 조직 공법과 면역치료 연구, 고품질 제품과 적절한 문서 준비가 필수적입니다. Gibco™ Cell Therapy Systems(CTS™) 브랜드에는 배지, 시약, 성장인자, 효소, 선택 비드, 장치 등 세포 치료 연구 분야에 사용할 수 있는 광범위한 고품질 제품이 포함되어 있습니다. 이 제품들은 21 CFR Part 820 품질 시스템 규정에 따라 제조되었으며 ISO 13485와 ISO 9001 인증을 받았습니다.

CTS 브랜드가 제공하는 장점:

cGMP-compliant documentation

- 시험성적서, 원산지 증명서, 약물 마스터 파일을 포함한 추적 문서
- 임상시험용 신약(IND)을 조정하는 시간 절감
- 모든 제품을 IND(CMC 가능)의 일부로 사용 가능
- CTS 제품 표시 및 용도 증명서
- cGMP및 FDA 준수 제품 제형

Seamless transition from research to clinic

- 합성 제형이므로 로드 간 편차최소화
- 확장 가능한 cGMP 조정에서 제조
- 상보성 RUO 제품
- 대부분의 제품에 멸균, 엔도독신, 외래 물질, 마이코플라즈마에 대한 포괄적인 QA 검사 실시

Expert consultation


- 모든 CTS 제품에 지역별 기술 지원 제공
- 노련한 전 세계 전문가가 연구부터 상용 단계까지 규제 관련 진행을 지원
- 모든 질문에 답변할 수 있는 세포치료 전문성
- 50년 넘게 축적된 Gibco™ 배지 경험

행사리스트

	Cat. no	Description	Size
Basal Media	11095080	MEM	500 mL
	11095088	MEM	10 X 500 mL
	12561056	MEM α, no Nucleosides	500 mL
	12561072	MEM α, no Nucleosides	10 X 500 mL
	12571063	MEM α, Nucleosides	500 mL
	12571071	MEM α, Nucleosides	10 X 500 mL
	11320033	DMEM/F-12	500 mL
	11320082	DMEM/F-12	10 X 500 mL
	11330032	DMEM/F-12, HEPES	500 mL
	11330057	DMEM/F-12, HEPES	10 X 500 mL
	11875093	RPMI 1640 Medium	500 mL
	11875119	RPMI 1640 Medium	10 X 500 mL
	22400089	RPMI 1640 Medium, HEPES	500 mL
	22400105	RPMI 1640 Medium, HEPES	10 X 500 mL
	11885084	DMEM, Low Glucose, Pyruvate	500 mL
	11885092	DMEM, Low Glucose, Pyruvate	10 X 500 mL
Balanced Salt Solution	11965092	DMEM, High Glucose	500 mL
	11965118	DMEM, High Glucose	10 X 500 mL
	11995065	DMEM, High Glucose, Pyruvate	500 mL
	11995073	DMEM, High Glucose, Pyruvate	10 X 500 mL
	12440053	Iscove S Modified Dulbecco Medium	500 mL
	12440061	Iscove S Modified Dulbecco Medium	10 X 500 mL
	10010023	PBS, pH 7.4	500 mL
	10010049	PBS, pH 7.4	10 X 500 mL
	20012027	PBS, pH 7.2	500 mL
	20012050	PBS, pH 7.2	10 X 500 mL
	14190144	DPBS, no calcium, no magnesium	500 mL
	14190250	DPBS, no calcium, no magnesium	10 X 500 mL
	14175095	HBSS, no calcium, no magnesium, no phenol red	500 mL
	14175103	HBSS, no calcium, no Magnesium, no Phenol Red	10 X 500 mL
	12604013	TrypLE™ Express (1X), no Phenol Red	100 mL
Dissociation Reagent	12604021	TrypLE™ Express (1X), no Phenol Red	500 mL
	12605010	TrypLE™ Express (1X), Phenol Red	100 mL
	12605036	TrypLE™ Express (1X), Phenol Red	20 X 100 mL
	15050065	0.25% Trypsin (1X), Phenol Red	100 mL
	15400054	0.5% Trypsin-EDTA (10X), no Phenol Red	10 X 100 mL
	25200056	0.25% Trypsin-EDTA (1X), Phenol Red	100 mL
	25200072	0.25% Trypsin-EDTA (1X), Phenol Red	500 mL
	25200114	0.25% Trypsin-EDTA (1X), Phenol Red	20 X 100 mL
	25300054	0.05% Trypsin-EDTA (1X), Phenol Red	100 mL
	25300062	0.05% Trypsin-EDTA (1X), Phenol Red	500 mL
	25300120	0.05% Trypsin-EDTA (1X), Phenol Red	20 X 100 mL
	15140122	Penicillin-Streptomycin (10,000 U/mL)	100 mL
	15140163	Penicillin-Streptomycin, Liquid (10,000 units Penicillin/10,000 µg Streptomycin)	20 X 100 mL
	15240062	Antibiotic-Antimycotic, 100X	100 mL
	15240112	Antibiotic-Antimycotic, 100X	20 X 100 mL
Antibiotics/ Antimycotics	15290018	Amphotericin B	20 mL
	15290026	Amphotericin B	50 mL
	15710072	Gentamicin (10 mg/ml)	10 X 10 mL
	10131027	Geneticin® Selective Antibiotic (G418 Sulfate)	100 mL
	10131035	Geneticin® Selective Antibiotic (G418 Sulfate)	20 mL
	11811023	Geneticin® Selective Antibiotic (G418 Sulfate), Powder	1 g
	11811031	Geneticin® Selective Antibiotic (G418 Sulfate), Powder	5g
	A1113902	Blasticidin S HCl (10 mg/mL)	10 X 1 mL
	A1113903	Blasticidin S HCl (10 mg/mL)	20 mL
	R21001	Blasticidin S HCl, Powder	50 mg
	R25001	Zeocin™ Selection Antibiotic (100 mg/mL)	1 g
	R25005	Zeocin™ Selection Antibiotic (100 mg/mL)	5 g
	A1113802	Puromycin Dihydrochloride	20 mL
	A1113803	Puromycin Dihydrochloride	10 X 1 mL
	10687010	Hygromycin B (50 mg/mL)	20 mL
Selection Antibiotics	15250061	Trypan Blue Solution, 0.4%	100 mL
	25030081	L-Glutamine, 200 mM Solution	100 mL
	25030164	L-Glutamine, 200 mM Solution	20 X 100 mL
	35050061	GlutaMAX™ Supplement	100 mL
	35050079	GlutaMAX™ Supplement	20 X 100 mL
	12648010	Recovery™ Cell Culture Freezing Medium	50 mL



	Cat. no	Description	Size
Gibco Recombinant Protein	13256029	hFGF Recombinant Human Protein	10 ug
	13257019	NGF Mouse Protein, Native, 2.5S Subunit	10 ug
	PHG0323	TNF-α Recombinant Human	1000 ug
	PHG3013	TNF-α Recombinant Human	1000 ug
	PHG3016	TNF-α Recombinant Human	50 ug
	PHG3991	VEGF Recombinant Human Protein	100 ug
	PHG3994	VEGF Recombinant Human Protein	10 ug
	PHG0021	FGF-Basic (AA 10-155) Recombinant Human Protein	100 ug
	PHG0023	FGF-Basic (AA 10-155) Recombinant Human Protein	1 mg
	PHG0024	FGF-Basic (AA 10-155) Recombinant Human Protein	10 ug
	PHG0026	FGF-Basic (AA 10-155) Recombinant Human Protein	50 ug
	PHG0041	PDGF-BB Recombinant Human Protein	100 ug
	PHG0043	PDGF-BB Recombinant Human Protein	1000 ug
	PHG0154	FGF4 Recombinant Human Protein	25 ug
	PHG0254	HGF Recombinant Human Protein	10 ug
	PHG0261	FGF-Basic (AA 1-155) Recombinant Human Protein	100 ug
	PHG0263	FGF-Basic (AA 1-155) Recombinant Human Protein	1 mg
	PHG0311	EGF Recombinant Human Protein	100 ug
	PHG0313	EGF Recombinant Human Protein	1 mg
	PHG0314	EGF Recombinant Human Protein	10 ug
	PHG9204	TGFB1 Recombinant Human Protein	5 ug
	PHG9214	TGFB1 Recombinant Human Protein	10 ug
	PMC2011	GM-CSF Recombinant Mouse Protein	100 ug
	PMC3016	TNF Recombinant Mouse Protein	25 ug
	PMG0034	FGF-Basic Recombinant Mouse Protein	10 ug
Transfection Reagent	03112	β Amyloid [1-42] PTD Human Protein	500 ug
	18324012	Lipofectamine® Transfection Reagent	1 mL
	18324020	Lipofectamine® Transfection Reagent	4 X 1 mL
	11668027	Lipofectamine® 2000 Transfection Reagent	0.75 mL
	11668019	Lipofectamine® 2000 Transfection Reagent	1.5 mL
	15338100	Lipofectamine® LTX & Plus Reagent	1 mL
	15338500	Lipofectamine® LTX & Plus Reagent	15 mL
	L3000008	Lipofectamine® 3000 Transfection Reagent	0.75 mL
	L3000015	Lipofectamine® 3000 Transfection Reagent	1.5 mL
	LMRNA003	Lipofectamine® MessengerMAX	0.3 mL
	LMRNA008	Lipofectamine® MessengerMAX	0.75 mL
	LMRNA015	Lipofectamine® MessengerMAX	1.5 mL
	11514015	PLUS™ Reagent	0.85 mL
	12252011	Oligofectamine™ Transfection Reagent	1 mL
	13778075	Lipofectamine® RNAiMAX Transfection Reagent	0.75 mL
	13778150	Lipofectamine® RNAiMAX Transfection Reagent	1.5 mL
Gateway Cloning	10362100	Cellfectin® II Reagent	1 mL
	16447100	FreeStyle™ MAX Reagent	1 mL
	A14524	ExpFectamine™ 293 Transfection Kit	1 kit
	11789013	Gateway® BP Clonase® Enzyme Mix	20 reactions
	11789020	Gateway® BP Clonase® II Enzyme mix	20 reactions
	11789100	Gateway® BP Clonase® II Enzyme mix	100 reactions
	11791019	Gateway® LR Clonase® Enzyme mix	20 reactions
	11791043	Gateway® LR Clonase® Enzyme Mix	100 reactions
	11791020	Gateway® LR Clonase® II Enzyme mix	20 reactions
	11791100	Gateway® LR Clonase® II Enzyme mix	100 reactions
Sample Prep	K240020	pENTR™/D-TOPO® Cloning Kit with One Shot® TOP10 Chemically Competent E. coli	20 reactions
	K240020SP	pENTR™/D-TOPO® Cloning Kit	20 reactions
	K2402020	pENTR™/SD/D-TOPO® Cloning Kit with One Shot® TOP10 Chemically Competent E. coli	20 reactions
	K250020	pCR4®/GW/TOPO® TA Cloning Kit	20 reactions
	K252020	pCR4®/B/GW/TOPO® TA Cloning Kit with One Shot® Mach1™/T1R Chemically Competent E. coli	20 reactions
	K210012	PureLink® Quick Gel Extraction Kit	50 preps
	K210025	PureLink® Quick Gel Extraction Kit	250 preps
	K220001	PureLink® Quick Gel Extraction and PCR Purification Combo Kit	50preps
	K310001	PureLink® PCR Purification Kit	50preps
	10603027	DNAzo® Reagent	100 mL
	K182001	PureLink® Genomic DNA Mini Kit	50 preps
	K182002	PureLink® Genomic DNA Mini Kit	250 preps
	K210001	PureLink® HQ Mini Plasmid DNA Purification Kit	100 preps
	K210002	PureLink® HiPure Plasmid Miniprep Kit	25 preps
	K210003	PureLink® HiPure Plasmid Miniprep Kit	100 preps
	K210004	PureLink® HiPure Plasmid Midiprep Kit	25 preps



	Cat. no	Description	Size
Sample Prep	K210010	PureLink® Quick Plasmid Miniprep Kit	50 preps
	K210011	PureLink® Quick Plasmid Miniprep Kit	250 preps
	K210014	PureLink® HiPure Plasmid Filter Midprep Kit	25 preps
	K210015	PureLink® HiPure Plasmid Filter Midprep Kit	50 preps
	K210016	PureLink® HiPure Plasmid Filter Midprep Kit	10 preps
	K210017	PureLink® HiPure Plasmid Filter Maxiprep Kit	25 preps
	K210018	PureLink® HiPure BAC buffer kit	1 kit
	AM1560	miVana™ miRNA Isolation Kit, with phenol	40 preps
	AM1561	miVana™ miRNA Isolation Kit, without phenol	40 preps
	AM1975	RecoverAll™ Total Nucleic Acid Isolation Kit for FFPE	40 preps
	10296010	TRIzol® LS Reagent	100 mL
	10296028	TRIzol® LS Reagent	200 mL
	12183555	TRIzol® Plus RNA Purification System	50 preps
	15596018	TRIzol® Reagent	200 mL
	15596026	TRIzol® Reagent	100 mL
	AM9738	TRI Reagent®	100 mL
	AM7020	RNAlater® RNA Stabilization Solution without manual	100 mL
	AM7021	RNAlater® RNA Stabilization Solution	500 mL
	AM7024	RNAlater® RNA Stabilization Solution	250 mL
	A13282	GeneAid® Site-Directed Mutagenesis System	16 reactions
	Q33217	Qubit® 3.0 Quantitation Starter Kit	1 kit
	Q33218	Qubit® 3.0 NGS Starter Kit	1 kit
	Q32850	Qubit® dsDNA BR Assay Kit	100 assays
	Q32851	Qubit® dsDNA HS Assay Kit	100 assays
	Q32852	Qubit® RNA HS Assay Kit	100 assays
	Q32853	Qubit® dsDNA BR Assay Kit	500 assays
	Q32854	Qubit® dsDNA HS Assay Kit	500 assays
	Q32855	Qubit® RNA HS Assay Kit	500 assays
	Q10210	Qubit® RNA BR Assay Kit	100 assays
	Q10211	Qubit® RNA BR Assay Kit	500 assays
	Q10212	Qubit® ssDNA Assay Kit	100 assays
	Q32880	Qubit® microRNA Assay Kit	100 assays
	Q32881	Qubit® microRNA Assay Kit	500 assays
	Q33211	Qubit® Protein Assay Kit	100 assays
	Q33212	Qubit® Protein Assay Kit	500 assays
	Q32856	Qubit® Assay Tubes	500 tubes
Antifades & Mountants	I36933	Image-It™ Fx Signal Enhancer	10ml
	P36930	Prolong Gold antifade reagent	10ml
	P36931	Prolong® Gold antifade reagent with DAPI	10ml
	P36934	Prolong® Gold antifade reagent - special packaging	5x2ml
	P36970	Prolong® Diamond Antifade Mountant	10 mL
	P36971	Prolong® Diamond Antifade Mountant with DAPI	10 mL
	P36961	Prolong® Diamond Antifade Mountant	5 x 2 mL
	P36962	Prolong® Diamond Antifade Mountant with DAPI	5 x 2 mL
	P36974	Prolong® Live Antifade Reagent	5 x 1 mL
	P36975	Prolong® Live Antifade Reagent	1 mL
	S36936	SlowFade® Gold antifade reagent	10ml
	S36938	SlowFade® Gold antifade reagent with DAPI	10ml
	S36937	SlowFade® Gold antifade reagent "special packaging"	5x2ml
	S36973	SlowFade® Diamond Antifade Mountant with DAPI	10 mL
	S36964	SlowFade® Diamond Antifade Mountant with DAPI	5 x 2 mL
Prism Reagent & Cell Viability Assay	R37605	NucBlue™ Live ReadyProbes™ Reagent	1 kit
	R37606	NucBlue® Fixed Cell ReadyProbes™ Reagent	1 kit
	R37602	Image-It™ Fixation/Permeabilization Kit	1 kit
	R37603	BackDrop® Background Suppressor ReadyProbes™ Reagent	1 kit
	A14291DJ	Live Cell Imaging Solution	500 ml
	A13261	PrestoBlue™ Cell Viability Reagent	25ml
	L3224	LIVE/DEAD® Viability/Cytotoxicity Kit, for mammalian cells	1 kit
	A22066	ATP Determination Kit	1 kit
	A11001	Alexa Fluor® 488 goat anti-mouse IgG (H+L) "2 mg/mL"	0.5ml
	A11005	Alexa Fluor® 594 goat anti-mouse IgG (H+L) "2 mg/mL"	0.5ml
	A21235	Alexa Fluor® 647 goat anti-mouse IgG (H+L) "2 mg/mL"	0.5ml
	A11009	Alexa Fluor® 488 goat anti-rabbit IgG (H+L) "2 mg/mL"	0.5ml
	A21428	Alexa Fluor® 555 goat anti-rabbit IgG (H+L) "2 mg/mL"	0.5ml
	A11012	Alexa Fluor® 594 goat anti-rabbit IgG (H+L) "2 mg/mL"	0.5ml
	A21244	Alexa Fluor® 647 goat anti-rabbit IgG (H+L) "2 mg/mL"	0.5ml
Secondary Antibody	A21422	Alexa Fluor® 555 goat anti-mouse IgG (H+L) "2 mg/mL"	0.5ml
	G21040	GOAT ANTI-MOUSE IGG (H+L) 1 MG	1 mg
	G21234	GOAT ANTI-RABBIT IGG (H+L) 1 MG	1 mg
Nupage Precasting Gel	EA03552BOX	NUPAGE® Novex® 7% Tris-Acetate Protein Gels, 1.0 mm, 12 well	10 ea/ box
	EA03554BOX	NUPAGE® Novex® 7% Tris-Acetate Protein Gels, 1.0 mm, 15 well	10 ea/ box
	EA0355BOX	NUPAGE® Novex® 7% Tris-Acetate Protein Gels, 1.0 mm, 10 well	10 ea/ box

	Cat. no	Description	Size
Nupage Precasting Gel	EA03752BOX	NUPAGE® Novex 3-8% Tris-Acetate Gel 1.0 mm, 12 Well	10 ea/ box
	EA03755BOX	NUPAGE® Novex 3-8% Tris-Acetate Gel 1.0 mm, 15 Well	10 ea/ box
	EA03759BOX	NUPAGE® Novex 3-8% Tris-Acetate Gel 1.0 mm, 10 Well	10 ea/ box
	NP03019BOX	NUPAGE® Novex 10% Bis-Tris Gels, 1.0 mm 10 well	10 ea/ box
	NP03219BOX	NUPAGE® Novex 4-12% Bis-Tris Gels, 1.0 mm 10 well	10 ea/ box
	NP03229BOX	NUPAGE® Novex 4-12% Bis-Tris Gels, 1.0 mm 12 well	10 ea/ box
	NP03239BOX	NUPAGE® Novex 4-12% Bis-Tris Gels, 1.0 mm 15 well	10 ea/ box
	NP03369BOX	NUPAGE® Novex 4-12% Bis-Tris Gels, 1.5 mm 10 well	10 ea/ box
	NP03369BOX	NUPAGE® Novex 4-12% Bis-Tris Gels, 1.5 mm 15 well	10 ea/ box
	NP03419BOX	NUPAGE® Novex 12% Bis-Tris Gels, 1.0 mm 10 well	10 ea/ box
	NP03429BOX	NUPAGE® Novex 12% Bis-Tris Gels, 1.0 mm 12 well	10 ea/ box
iBlot Transfer Stack	IB301001	iBlot® Transfer Stack, NC Regular	10 ea
	IB301002	iBlot® Transfer Stack, NC Mini	10 ea
	IB301031	iBlot® Transfer Stack, NC Regular	3 X 10 ea
	IB301032	iBlot® Transfer Stack, NC Mini	3 X 10 ea
	IB401001	iBlot® Transfer stack, PVDF Regular	10 ea
	IB401002	iBlot® Transfer stack, PVDF Mini	10 ea
	IB401031	iBlot® Transfer stack, PVDF Regular	3 X 10 ea
	IB401032	iBlot® Transfer stack, PVDF Mini	3 X 10 ea
Protein Electrophoresis	E10001	XCell SureLock™ Mini-Cell	1 ea
	E10002	XCell SureLock™ Mini-Cell and XCell II™ Blot Module	1 kit
	LC2675	Novex® Tris-Glycine SDS Running Buffer (10X)	500 mL
	LC2676	Novex® Tris-Glycine SDS Sample Buffer (2X)	200 mL
	LC3675	Novex® Tris-Glycine Transfer Buffer (25X)	500 mL
	NP0001	NUPAGE® Mops SDS Running Buffer	500 mL
	NP0002	NUPAGE® MES SDS Running Buffer	500 mL
	NP0004	NUPAGE® sample reducing Buffer	250 uL
	NP00061	NUPAGE® Transfer Buffer (20X)	1 L
	B0001	Bolt® MOPS SDS Running Buffer (20X)	500 mL
	B0002	Bolt® MES SDS Running Buffer (20X)	500 mL
	B0004	10X Bolt® Sample Reducing Agent	250 uL
	B0007	4X Bolt™ LDS Sample Buffer	10 mL
	BT00061	Bolt® Transfer Buffer (20X)	1L
	NP0007	NUPAGE® LDS Sample Buffer (4X)	10 mL
ELISA	LC5602	MagicMark™ XP Western Protein Standard	250 uL
	LC5800	Novex® Sharp Pre-stained Protein Standards	2 X 250 uL
	LC5925	SeeBlue® Plus2 Pre-Stained Standard	500 uL
	KH-B3442	Aβ42 Human ELISA Kit	192 tests
	KH-B3482	Aβ40 Human ELISA Kit	192 tests
	KH-B3544	Aβ 42 Human Ultrasensitive ELISA Kit	96 tests
	KH-C3011	TNF alpha ELISA Kit, Human	96 tests
	KMB3441	Aβ 42 Mouse ELISA Kit	96 tests
	KMC3011	TNF-α Mouse ELISA Kit	96 tests
	002023	TMB Chromogen Solution (for ELISA)	500 mL
	FNNO011	Cell Extraction Buffer	100 mL
	FNNO021	NP40 Cell Lysis Buffer	100 mL
Reverse Transcriptase	18057018	M-MLV Reverse Transcriptase Buffer	1 mL
	28025013	M-MLV Reverse Transcriptase	40,000 units
	28025021	M-MLV Reverse Transcriptase	200,000 units
	18060051	SuperScript® III First-Strand Synthesis System	50 reactions
	11904018	SuperScript® First-Strand Synthesis System for RT-PCR	50 reactions
	18064014	SuperScript® II Reverse Transcriptase	10,000 units
	18064022	SuperScript® II Reverse Transcriptase	2,000 units
	18064071	SuperScript® II Reverse Transcriptase	4 x 10,000 units
	18060044	SuperScript® III Reverse Transcriptase	10,000 units
	18060085	SuperScript® III Reverse Transcriptase	4 x 10,000 units
	18060093	SuperScript® III Reverse Transcriptase	2,000 units
	11754050	SuperScript® VLO™ cDNA Synthesis Kit	50 x 20 µL reactions
	11754250	SSuperScript® VLO™ cDNA Synthesis Kit	250 x 20 µL reactions
	18090050	Superscript® IV Reverse Transcriptase	10,000 units
	18090020	Superscript® IV Reverse Transcriptase	4 x 10,000 units
Cell Culture Ware	18091050	Superscript® IV First-Strand Synthesis System	50 reactions
	18091200	Superscript® IV First-Strand Synthesis System	200 reactions
	A70551	Nunc Cell Culture Treated EasyFlasks T25 Filter cap, #156367	200/cs
	A70555	Nunc Cell Culture Treated EasyFlasks T75 Filter cap, #156499	100/cs
	A70571	Nunc Cell Culture Treated EasyFlasks T175 Filter cap, #156910	30/cs
	A70550	Nunc Cell Culture Treated EasyFlasks T25 Vent/close cap, #156340	200/cs
	A70553	Nunc Cell Culture Treated EasyFlasks T75 Vent/close cap, #156472	100/cs
	A70573	Nunc Cell Culture Treated EasyFlasks T175 Vent/close cap, #156920	30/cs
	A70557	Nunc Cell Culture Treated EasyFlasks T175 Vent/Close cap, #156602	32/cs
	A70692	Nunc Cell Culture Treated Straight Flasks T175 Filter cap, #178883	32/cs
	A70395	Nunc Cell Culture Treated TripleFlasks Filter cap, #132913	32/cs



	Cat. no	Description	Size
Cell Culture Ware	A70516	Nunc Cell Culture/Petri Dishes 35mm, #153066	500/cs
	A70489	Nunc Cell Culture/Petri Dishes 60mm, #150289	400/cs
	A27721	Nunc Cell Culture/Petri Dishes 100mm, #172931	240/cs
	A70614	Nunc Cell Culture/Petri Dishes 150mm, #168381	80/cs
	A70676	Nunc Cell-Culture Treated Multidishes 4well, #176740	120/cs
	A70439	Nunc Cell-Culture Treated Multidishes 6well, #140675	75/cs
	A70497	Nunc Cell-Culture Treated Multidishes 12well, #150628	75/cs
	A70451	Nunc Cell-Culture Treated Multidishes 24well, #142475	75/cs
	A70499	Nunc Cell-Culture Treated Multidishes 48well, #150687	75/cs
	A70604	Nunc MicroWell 96-Well Microplates Cell-culture Treated, Flat bottom, #167008	50/cs
	A70582	Nunc MicroWell 96-Well Microplates Cell-culture Treated, Round bottom, #163320	50/cs
	A70404	Nunc F96 MicroWell Black Polystyrene Plate Cell-culture Treated, #137101	50/cs
	A70399	Nunc F96 MicroWell White Polystyrene Plate Cell-culture Treated, #136101	50/cs
	A67919	Nalgene Rapid-Flow Bottle Top Filters, PES membrane, 500mL, 45mm neck size, #595-4520	12/cs
Immunoassay	A67918	Nalgene Rapid-Flow Bottle Top Filters, PES membrane, 500mL, 35mm neck size, #595-3320	12/cs
	A67924	Nalgene Rapid-Flow Bottle Top Filters, PES membrane, 1000mL, 45mm neck size, #597-4520	12/cs
	A67920	Nalgene Rapid-Flow Filter Units, PES Membrane, 500mL, 0.2um pore size, #566-0020	12/cs
	A67921	Nalgene Rapid-Flow Filter Units, PES Membrane, 1000mL, 0.2um pore size, #567-0020	12/cs
	A68010	Nalgene Rapid-Flow Filter Units, PES Membrane, 500mL, 0.2um pore size, #569-0020	12/cs
	A67970	NalgeneRapid-FlowSterile Filter Storage Bottles, 250mL, #455-0250	24/cs
	A67923	NalgeneRapid-FlowSterile Filter Storage Bottles, 500mL, #455-0500	12/cs
	A67922	NalgeneRapid-FlowSterile Filter Storage Bottles, 1000mL, #455-1000	12/cs
	A71121	Nunc 96 Immunoplate, Maxiisorp, Certified, #439454	60/cs
	A71049	NuncCell Culture Cryogenic Tubes, 1.8mL external, Self-standing, #375418	450/cs
	A71000	NuncCell Culture Cryogenic Tubes, 1.8mL Internal, Self-standing, #368632	450/cs
	A71060	NuncCell Culture Cryogenic Tubes, 1.8mL Internal, Self-standing, starfoot, #377627	450/cs
	A71056	NuncTube/Vial Containers 4*10, #376589	1/cs
	A68060	NalgeneGeneral Long-Term Storage Cryogenic Tubes, 2.0mL external, #5000-0020	500/cs
Cryoware	A68979	Nalgene CryoCane, #5015-0001	12/pk
	A67098	Nalgene Polycarbonate CryoBox, 9*9 , #5026-0909	4/pk
	A68897	Thermo ScientificMr. FrostyFreezing Container for 3.6mL tubes, #5100-0036	1/unit
	A68898	Thermo ScientificMr. FrostyFreezing Container for 4.5mL tubes, #5100-0050	1/unit
	A68969	Thermo ScientificMr. FrostyFreezing Container for 1.0-2.0mL tubes, #5100-0001	1/unit
	A68910	Thermo ScientificBenchtop Coolers -20°C, 4*8, #5115-0032	1/unit
	A70697	NuncCell Scrapers, 23cm, #179693	50/cs
	A70699	NuncCell Scrapers, 32cm, #179707	50/cs
	A70740	NuncSealing Tapes, Cell culture, sterile applications, #236366	200/cs
	A70789	NuncDisposable Loops, 1*1L, #253287	1000/cs
	A67600	NalgeneWide-Mouth Straight-Sided PPCC Jars with Closure, 1000mL , #2118-0032	6/pk
	A67729	NalgeneAutoactivable Polypropylene Carboy with Spigot, 20L , #2319-0050	1/pk
	A68734	NalgeneOak Ridge High-Speed PPCC Centrifuge Tubes, 50mL , #3119-0050	10/pk
	A67311	NalgenePPCC Centrifuge Bottles, 250mL , #3120-0250	4/pk
Labware	A68041	NalgenePolypropylene Quick Action Spigots for Nalgene carboy, #6422-0010	1/pk
	A69623	QSP 1-200uL Gel loading tip, bulk, #1010-Q	1000/pk
	A69624	QSP 1-200uL Gel loading tip, racked, #1010-R204-Q	1020/pk
	A69256	QSP 10uL Micro tip, bulk, #104-Q	1000/pk
	A69561	QSP 10uL tip, reload #T104RL-Q	960/pk
	A69288	QSP 200uL Yellow tip, bulk, #110-Q	1000/pk
	A69334	QSP 200uL Yellow tip, graduated, bulk, #113-G-Q	1000/pk
	A69606	QSP 200uL Clear tip, thin wall, bulk , #TW110-N-Q	1000/pk
	A69569	QSP 200uL Yellow tip, beveled, reload, #T110BRL-Q	960/pk
	A69585	QSP 200uL Clear tip, beveled, graduated, reload, #T113RLN-Q	960/pk
	A69592	QSP 200uL Clear tip, thin wall, clear, reload #TTW110RLN-Q	960/pk
	A69305	QSP 1000uL Blue tip, bulk, #111-Q	1000/pk
	A69316	QSP 1250uL Clear tip, extra long, bulk, #112NXL-Q	1000/pk
	A69323	QSP 1250uL Blue tip, extra long, bulk, #112XL-Q	1000/pk
Tip & Tube	A69577	QSP 1250 uL Clear tip, extra long, graduated, reload, #T112NXLRL-Q	960/pk
	A69579	QSP 1250uL Clear tip, extra long, racked, #T112NXLRL-Q	960/pk
	A69586	QSP 10uL Filter tip, sterile, clear, #TF104-10-Q	960/pk
	A69634	QSP 2-20uL Filter tip, sterile, clear, #TF113-20-Q	960/pk
	A69590	QSP 200uL Filter tip, sterile, clear, #TF140-200-Q	960/pk
	A69587	QSP 1000uL Filter tip, sterile, clear, #TF112-1000-Q	960/pk
	A69091	QSP 0.2mL Thin wall tube with flat cap, #431-Q	1000/pk
	A69159	QSP 1.5mL Graduated microcentrifuge tubes, #509-GRD-Q	500/pk
	A69927	ART 10uL Filtered tip, #2139-HR	960/pk
	A69934	ART 10uL Filtered tip, extended length, #2140-HR	960/pk
	A69951	ART 20uL Filtered tip, #2149-HR	960/pk
	A69904	ART 200uL Filtered tip, #2069-HR	960/pk
	A69963	ART 1000uL Filtered tip, #2179-HR	768/pk
	A20480	ART 1000uL Filtered tip, extended length, #2079-HR	768/pk
A70341	ART 10uL, tp, bulk, #3500	1000/pk	
A19740	ART 10uL tip, non-sterile, racked, #3502-HR	960/pk	