


ThermoFisher

ThermoFisher 장비 구매 찬스

가격 할인 및 사은증정 행사



행사기간 : ~ 2013년 1월 31일까지

CureBio

(주) 큐어바이오

T. 02-427-3669 www.curebio.co.kr

장비도 구입하고 아이패드도 내손에

쏘~옥



5백만원 미만 : 가격 할인

5백만원 이상 : 가격 할인 및 아이패드미니 32G (또는 아이패드 4, 16G)

천만원 이상 : 가격 할인 및 아이패드미니 32G 2대 (또는 아이패드 4, 16G 2대)

천오백만원 이상 : 가격 할인 및 아이패드미니 32G 3대
(또는 아이패드 4, 16G, 3대)

** 사은품은 결제시 지급



Microcentrifuge & Megafuge

Cat. No.	Model	Max. RCF	Electrical	냉장유무
75002410	Heraeus Pico 17 microcentrifuge, includes 24 x 1.5/2.0 mL rotor	17,000 xg	230V 50/60Hz	비냉장
75002420	Heraeus Fresco 17 microcentrifuge, includes 24 x 1.5/2.0 mL rotor	17,000 xg	230V 50/60Hz	냉장
75002415	Heraeus Pico 21 microcentrifuge, includes 24 x 1.5/2.0 mL rotor	21,100 xg	230V 50/60Hz	비냉장
75002425	Heraeus Fresco 21 microcentrifuge, includes 24 x 1.5/2.0 mL rotor	21,100 xg	230V 50/60Hz	냉장

Cat. No.	Model	Swinging bucket rotors		Fixed angle rotors		냉장유무
		Max capacity	Max speed	Max capacity	Max speed	
75004230	Heraeus Megafuge 16	4 X 400ml	5,500 rpm	6 x 100 ml	15,200 rpm	비냉장
75004270	Heraeus Megafuge 16R	4 X 400ml	5,500 rpm	6 x 100 ml	15,200 rpm	냉장
75004503	Heraeus Megafuge 40	4 X 750ml	5,300 rpm	6 x 100 ml	15,200 rpm	비냉장
75004518	Heraeus Megafuge 40R	4 X 750ml	5,300 rpm	6 x 100 ml	15,200 rpm	냉장



CO2 Incubator

Cat. No.	Model	Interior	Volume
51023121	Heraeus BB15	stainless steel	151 L
51026280	Heracell 150i single chamber	stainless steel	150 L
51026281		100% pure copper	
51026333	Heracell 240i single chamber	stainless steel	240 L
51026334		100% pure copper	



Cryogenic System

Cat. No.	Model	Capacity	2ml Vial rack	box/rack	vial/box	Vial capacity 2ml
CY50945-70	Locator 8 Plus	121 L	8	10	25	2,000
CY50925-70	Locator Jr Plus	71 L	4	5	100	2,000
CY50935-70	Locator 4 Plus	121 L	4	10	100	4,000
CY50985-70	Locator 6 Plus	184 L	6	10	100	6,000



Ultra Low Temperature Upright Freezer

Cat. No.	Model	Capacity	Temp. range	Electrical
IE2386D	Isotemp basic ULT upright freezer	651.3 L	-50°C to -86°C	208-230V 60Hz
IE1786D	Isotemp basic ULT upright freezer	489.9 L	-50°C to -86°C	208-230V 60Hz
IU2886D	Isotemp ULT upright freezer	793 L	-50°C to -86°C	208-230V 61Hz



Biological Safety Cabinet

Cat. No.	Type	Exterior Dimensions H x W x D (mm)	Interior Dimensions H x W x D (mm)
1376	1300 Series Class II, Type A2	1568 x 1300 x 800	780 x 1200 x 465
1378	1300 Series Class II, Type A2	1568 x 1900 x 800	780 x 1800 x 465